

产品名称: Nestin Rabbit Monoclonal Antibody
产品货号: AMRe87738



产品概述 (Summary)

产品名称 (Production Name)	Nestin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Nestin
别名 (Alternative Names)	RC2; Marc2; ESTM46; Ifaprc2
基因 ID (Gene ID)	18008
蛋白 ID (SwissProt ID)	Q6P5H2, P21263.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW;; Observed MW:

研究背景 (Background)


产品名称: Nestin Rabbit Monoclonal Antibody
产品货号: AMRe87738



Predicted to enable CCR5 chemokine receptor binding activity. Acts upstream of or within cell projection morphogenesis; negative regulation of neuron apoptotic process; and nervous system development. Located in cytoplasm and intermediate filament. Is expressed in several structures, including central nervous system; embryo mesenchyme; eye; foregut-midgut junction; and metanephros. Orthologous to human NES (nestin). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Neuro-2a cells using AMRe87738 at 1:1000.

注意事项 (Note)

For research use only .