

产品名称: Cullin 3 Rabbit Monoclonal Antibody
产品货号: AMRe87740



产品概述 (Summary)

产品名称 (Production Name)	Cullin 3 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cullin 3
别名 (Alternative Names)	CUL-3; PHA2E
基因 ID (Gene ID)	8452
蛋白 ID (SwissProt ID)	Q13618.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:89 kDa; Observed MW:89 kDa

研究背景 (Background)


产品名称: Cullin 3 Rabbit Monoclonal Antibody
产品货号: AMRe87740




This gene encodes a member of the cullin protein family. The encoded protein plays a critical role in the polyubiquitination and subsequent degradation of specific protein substrates as the core component and scaffold protein of an E3 ubiquitin ligase complex. Complexes including the encoded protein may also play a role in late endosome maturation. Mutations in this gene are a cause of type 2E pseudohypoaldosteronism. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Mar 2012]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of Cullin 3/CUL-3 in K562, 3T3, HeLa cell lysates using Cullin 3/CUL-3 antibody (1:1000 diluted).



Western blot analysis of extracts from C2C12 cells using AMRe87740 at 1:5000.

注意事项 (Note)

For research use only .