

产品名称: GC1q R Rabbit Monoclonal Antibody
产品货号: AMRe87743



产品概述 (Summary)

产品名称 (Production Name)	GC1q R Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GC1q R
别名 (Alternative Names)	p32; HABP1; gC1qR; GC1QBP; SF2p32; gC1Q-R
基因 ID (Gene ID)	708
蛋白 ID (SwissProt ID)	Q07021.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:31 kDa; Observed MW:31 kDa

研究背景 (Background)


产品名称: GC1q R Rabbit Monoclonal Antibody
产品货号: AMRe87743




The human complement subcomponent C1q associates with C1r and C1s in order to yield the first component of the serum complement system. The protein encoded by this gene is known to bind to the globular heads of C1q molecules and inhibit C1 activation. This protein has also been identified as the p32 subunit of pre-mRNA splicing factor SF2, as well as a hyaluronic acid-binding protein. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of GC1q R in K562, Rat Brain, C6, 3T3, Hela cell lysates using GC1q R antibody (1:1000 diluted).



Western blot analysis of extracts from HCT116 cells and Rat stomach, Mouse kidney tissue using AMRe87743 at 1:4000.

注意事项 (Note)

For research use only .