

产品名称: CPSF6 Rabbit Monoclonal Antibody
产品货号: AMRe87746



产品概述 (Summary)

产品名称 (Production Name)	CPSF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CPSF6
别名 (Alternative Names)	CFIM; CFIM68; CFIM72; HPBR11-4; HPBR11-7
基因 ID (Gene ID)	11052
蛋白 ID (SwissProt ID)	Q16630.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:20000,IHC 1:100-1:200,ICC/IF 1:100-1:200,FC 1:10-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:59 kDa; Observed MW:70 kDa

产品名称: CPSF6 Rabbit Monoclonal Antibody
产品货号: AMRe87746




研究背景 (Background)


The protein encoded by this gene is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of the protein with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors. The cleavage factor complex is composed of four polypeptides. This gene encodes the 68kD subunit. It has a domain organization reminiscent of spliceosomal proteins. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot detection of CPSF6 in Rat Brain, MCF7, 3T3, C6 cell lysates using CPSF6 antibody (1:1000 diluted).



Western blot analysis of extracts from HepG2 cells and Mouse heart tissue using AMRe87746 at 1:5000.

注意事项 (Note)

For research use only .