

产品名称: ALKBH5 Rabbit Monoclonal Antibody
产品货号: AMRe87757



产品概述 (Summary)

产品名称 (Production Name)	ALKBH5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALKBH5
别名 (Alternative Names)	ABH5; OFOX D; OFOX D1
基因 ID (Gene ID)	54890
蛋白 ID (SwissProt ID)	Q6P6C2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:44 kDa; Observed MW:44-48 kDa

研究背景 (Background)


产品名称: ALKBH5 Rabbit Monoclonal Antibody
产品货号: AMRe87757




Dioxygenase that demethylates RNA by oxidative demethylation: specifically demethylates N6-methyladenosine (m6A) RNA, the most prevalent internal modification of messenger RNA (mRNA) in higher eukaryotes (PubMed:23177736, PubMed:24778178, PubMed:24616105, PubMed:24489119). Can also demethylate N6-methyladenosine in single-stranded DNA (in vitro) (PubMed:24616105). Requires molecular oxygen, alpha-ketoglutarate and iron (PubMed:21264265, PubMed:23177736, PubMed:24778178, PubMed:24616105, PubMed:24489119). Demethylation of m6A mRNA affects mRNA processing and export (PubMed:23177736). Required for spermatogenesis (By similarity).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87757 at 1:1000.



Western blot analysis of extracts from HepG2, C2C12, PC-12 cells using AMRe87757 at 1:1000.

注意事项 (Note)

For research use only .