

产品名称: CLN5 Rabbit Monoclonal Antibody
产品货号: AMRe87761



产品概述 (Summary)

产品名称 (Production Name)	CLN5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLN5
别名 (Alternative Names)	NCL
基因 ID (Gene ID)	1203
蛋白 ID (SwissProt ID)	O75503.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:2000,ICC/IF 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:42 kDa; Observed MW:50 kDa

研究背景 (Background)


产品名称: CLN5 Rabbit Monoclonal Antibody
产品货号: AMRe87761



This gene is one of eight which have been associated with neuronal ceroid lipofuscinoses (NCL). Also referred to as Batten disease, NCL comprises a class of autosomal recessive, neurodegenerative disorders affecting children. The genes responsible likely encode proteins involved in the degradation of post-translationally modified proteins in lysosomes. The primary defect in NCL disorders is thought to be associated with lysosomal storage function.[provided by RefSeq, Oct 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from SCC-9 cells using db11462 at 1:3000.

注意事项 (Note)

For research use only .