

产品名称: RICTOR Rabbit Monoclonal Antibody
产品货号: AMRe87778



产品概述 (Summary)

产品名称 (Production Name)	RICTOR Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,FC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RICTOR
别名 (Alternative Names)	PIA; AVO3; hAVO3
基因 ID (Gene ID)	253260
蛋白 ID (SwissProt ID)	Q6R327.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,FC 1:200-1:500
蛋白分子量 (Molecular Weight)	Calculated MW:192 kDa; Observed MW:200 kDa

研究背景 (Background)


产品名称: RICTOR Rabbit Monoclonal Antibody
产品货号: AMRe87778




RICTOR and MTOR (FRAP1; MIM 601231) are components of a protein complex that integrates nutrient- and growth factor-derived signals to regulate cell growth (Sarbassov et al., 2004 [PubMed 15268862]).[supplied by OMIM, Mar 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from HeLa cells using AMRe87778 at 1:1000.



Western blot analysis of extracts from SW620, Caco2, 293T cells using AMRe87778 at 1:1000.

注意事项 (Note)

For research use only .