

产品名称: TIMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87792



产品概述 (Summary)

产品名称 (Production Name)	TIMP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TIMP1
别名 (Alternative Names)	EPA; Clgi; Timp; TIMP-1; TPA-S1
基因 ID (Gene ID)	21857
蛋白 ID (SwissProt ID)	P12032.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW:23 kDa; Observed MW:26 kDa

研究背景 (Background)


产品名称: TIMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87792



Metalloproteinase inhibitor that functions by forming one to one complexes with target metalloproteinases, such as collagenases, and irreversibly inactivates them by binding to their catalytic zinc cofactor. Acts on MMP1, MMP2, MMP3, MMP7, MMP8, MMP9, MMP10, MMP11, MMP12, MMP13 and MMP16. Does not act on MMP14 (By similarity). Also functions as a growth factor that regulates cell differentiation, migration and cell death and activates cellular signaling cascades via CD63 and ITGB1. Plays a role in integrin signaling.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from C2C12 cells using AMRe87792 at 1:1000.

注意事项 (Note)

For research use only .