

产品名称: VPS24 Rabbit Monoclonal Antibody
产品货号: AMRe87794



产品概述 (Summary)

产品名称 (Production Name)	VPS24 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	VPS24
别名 (Alternative Names)	NEDF; VPS24; CGI-149
基因 ID (Gene ID)	51652
蛋白 ID (SwissProt ID)	Q9Y3E7.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:50-1:100,FC 1:10-1:100,IP 1:10-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:25 kDa; Observed MW:25 kDa

研究背景 (Background)


产品名称: VPS24 Rabbit Monoclonal Antibody
产品货号: AMRe87794



This gene encodes a protein that sorts transmembrane proteins into lysosomes/vacuoles via the multivesicular body (MVB) pathway. This protein, along with other soluble coiled-coil containing proteins, forms part of the ESCRT-III protein complex that binds to the endosomal membrane and recruits additional cofactors for protein sorting into the MVB. This protein may also co-immunoprecipitate with a member of the IFG-binding protein superfamily. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the upstream ring finger protein 103 (RNF103) gene. [provided by RefSeq, Nov 2010]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Caco2, C2C12 cells using AMRe87794 at 1:3000.

注意事项 (Note)

For research use only .