

产品名称: Gephyrin Rabbit Monoclonal Antibody
产品货号: AMRe87797



产品概述 (Summary)

产品名称 (Production Name)	Gephyrin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Gephyrin
别名 (Alternative Names)	GPH; GEPH; HKPX1; GPHRYN; MOCODC
基因 ID (Gene ID)	10243
蛋白 ID (SwissProt ID)	Q9NQX3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:80 kDa; Observed MW:93 kDa

研究背景 (Background)


产品名称: Gephyrin Rabbit Monoclonal Antibody
产品货号: AMRe87797




This gene encodes a neuronal assembly protein that anchors inhibitory neurotransmitter receptors to the postsynaptic cytoskeleton via high affinity binding to a receptor subunit domain and tubulin dimers. In nonneuronal tissues, the encoded protein is also required for molybdenum cofactor biosynthesis. Mutations in this gene may be associated with the neurological condition hyperplexia and also lead to molybdenum cofactor deficiency. Numerous alternatively spliced transcript variants encoding different isoforms have been described; however, the full-length nature of all transcript variants is not currently known. [provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from 293T cells using AMRe87797 at 1:1000.



Western blot analysis of extracts from HepG2, HCT116 cells and Rat liver tissue using AMRe87797 at 1:1000.

注意事项 (Note)

For research use only .