

产品名称: Cardiac Troponin I Rabbit Monoclonal Antibody
产品货号: AMRe87800



产品概述 (Summary)

产品名称 (Production Name)	Cardiac Troponin I Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Cardiac Troponin I
别名 (Alternative Names)	CMH7; RCM1; cTnI; CMD2A; TNNC1; CMD1FF
基因 ID (Gene ID)	7137
蛋白 ID (SwissProt ID)	P19429.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:50000,IHC 1:200-1:500,FC 1:20-1:50,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW:24 kDa; Observed MW:28 kDa

研究背景 (Background)


产品名称: Cardiac Troponin I Rabbit Monoclonal Antibody
产品货号: AMRe87800



Troponin I (TnI), along with troponin T (TnT) and troponin C (TnC), is one of 3 subunits that form the troponin complex of the thin filaments of striated muscle. TnI is the inhibitory subunit; blocking actin-myosin interactions and thereby mediating striated muscle relaxation. The TnI subfamily contains three genes: TnI-skeletal-fast-twitch, TnI-skeletal-slow-twitch, and TnI-cardiac. This gene encodes the TnI-cardiac protein and is exclusively expressed in cardiac muscle tissues. Mutations in this gene cause familial hypertrophic cardiomyopathy type 7 (CMH7) and familial restrictive cardiomyopathy (RCM).
[provided by RefSeq, Jul 2008]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Human heart tissue using AMRe87800 at 1:1000.

注意事项 (Note)

For research use only .