

产品名称: RAB13 Rabbit Monoclonal Antibody  
产品货号: AMRe87802



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | RAB13 Rabbit Monoclonal Antibody       |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IP                                  |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                                                                                              |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                                                                                 |
| 修饰 (Modification)   | Unmodified                                                                                                                                                   |
| 同种型 (Isotype)       | IgG                                                                                                                                                          |
| 克隆 (Clonality)      | Monoclonal                                                                                                                                                   |
| 形式 (Form)           | Liquid                                                                                                                                                       |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.                                                                     |
| 储存溶液 (Buffer)       | Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt. |
| 纯化方式 (Purification) | Affinity Purification                                                                                                                                        |

## 免疫原信息 (Immunogen)

|                        |         |
|------------------------|---------|
| 基因名 (Gene Name)        | RAB13   |
| 别名 (Alternative Names) | GIG4    |
| 基因 ID (Gene ID)        | 5872    |
| 蛋白 ID (SwissProt ID)   | P51153. |

## 产品应用 (Application)

|                          |                                          |
|--------------------------|------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IP 1:20-1:50             |
| 蛋白分子量 (Molecular Weight) | Calculated MW:23 kDa; Observed MW:23 kDa |

## 研究背景 (Background)


**产品名称: RAB13 Rabbit Monoclonal Antibody**  
**产品货号: AMRe87802**



This gene is a member of the Rab family of small G proteins and plays a role in regulating membrane trafficking between trans-Golgi network (TGN) and recycling endosomes (RE). The encoded protein is involved in the assembly of tight junctions, which are components of the apical junctional complex (AJC) of epithelial cells. The AJC plays a role in forming a barrier between luminal contents and the underlying tissue. Additional functions associated with the protein include endocytic recycling of occludin, regulation of epithelial cell scattering, neuronal regeneration and regulation of neurite outgrowth. Alternately spliced transcript variants have been observed for this gene. A pseudogene associated with this gene is located on chromosome 12. [provided by RefSeq, Jan 2013]

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of extracts from SW480 cells using AMRe87802 at 1:1000.

## 注意事项 (Note)

For research use only .