

产品名称: Vaspin Rabbit Monoclonal Antibody
产品货号: AMRe87823



产品概述 (Summary)

产品名称 (Production Name)	Vaspin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Vaspin
别名 (Alternative Names)	4632419J12Rik
基因 ID (Gene ID)	68054
蛋白 ID (SwissProt ID)	Q7TMF5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:48 kDa; Observed MW:48 kDa

研究背景 (Background)


产品名称: Vaspin Rabbit Monoclonal Antibody
产品货号: AMRe87823



Predicted to enable serine-type endopeptidase inhibitor activity. Acts upstream of or within negative regulation of gluconeogenesis; positive regulation of signal transduction; and regulation of lipid metabolic process. Located in extracellular space and plasma membrane. Is expressed in back skin. Orthologous to human SERPINA12 (serpin family A member 12). [provided by Alliance of Genome Resources, Apr 2022]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse stomach tissue using AMRe87823 at 1:1000.

注意事项 (Note)

For research use only .