

产品名称: FCGRT/FCRN Rabbit Monoclonal Antibody
产品货号: AMRe87827



产品概述 (Summary)

产品名称 (Production Name)	FCGRT/FCRN Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FCGRT/FCRN
别名 (Alternative Names)	FCRN; alpha-chain
基因 ID (Gene ID)	2217
蛋白 ID (SwissProt ID)	P55899.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:500-1:2000,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:40 kDa; Observed MW:40 kDa

研究背景 (Background)


产品名称: FCGRT/FCRN Rabbit Monoclonal Antibody
产品货号: AMRe87827



This gene encodes a receptor that binds the Fc region of monomeric immunoglobulin G. The encoded protein transfers immunoglobulin G antibodies from mother to fetus across the placenta. This protein also binds immunoglobulin G to protect the antibody from degradation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from THP-1 cells using AMRe87827 at 1:1000.

注意事项 (Note)

For research use only .