

产品名称: Folate Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe87835



产品概述 (Summary)

产品名称 (Production Name)	Folate Binding Protein Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Folate Binding Protein
别名 (Alternative Names)	FBP; FOLR
基因 ID (Gene ID)	2348
蛋白 ID (SwissProt ID)	P15328.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000
蛋白分子量 (Molecular Weight)	Calculated MW:30 kDa; Observed MW:40 kDa

产品名称: Folate Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe87835




研究背景 (Background)


The protein encoded by this gene is a member of the folate receptor family. Members of this gene family bind folic acid and its reduced derivatives, and transport 5-methyltetrahydrofolate into cells. This gene product is a secreted protein that either anchors to membranes via a glycosyl-phosphatidylinositol linkage or exists in a soluble form. Mutations in this gene have been associated with neurodegeneration due to cerebral folate transport deficiency. Due to the presence of two promoters, multiple transcription start sites, and alternative splicing, multiple transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Oct 2009]

研究领域 (Research Area)

图片 (Image Data)




Western blot detection of Folate Binding Protein in HeLa cell lysates using Folate Binding Protein antibody(1:1000 diluted).



Immunohistochemical analysis of paraffin-embedded human lung cancer using AMRe87835 antibody.

产品名称: Folate Binding Protein Rabbit Monoclonal Antibody
产品货号: AMRe87835



Western blot analysis of extracts from HeLa cells and Mouse kidney tissue using AMRe87835 at 1:1000.

注意事项 (Note)

For research use only .