

产品名称: LAMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87845



产品概述 (Summary)

产品名称 (Production Name)	LAMP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,FC,IP
种属反应性 (Reactivity)	Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% sodium azide and 0.05% protective protein. Stable for 12 months from date of receipt.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LAMP1
别名 (Alternative Names)	P2B; Perk; LGP-A; CD107a; Lamp-1; LGP-120
基因 ID (Gene ID)	16783
蛋白 ID (SwissProt ID)	P11438.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:200-1:500,ICC/IF 1:100-1:200,FC 1:200-1:500,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW:44 kDa; Observed MW:90-120 kDa

产品名称: LAMP1 Rabbit Monoclonal Antibody
产品货号: AMRe87845




研究背景 (Background)


Presents carbohydrate ligands to selectins. Also implicated in tumor cell metastasis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of extracts from Mouse liver tissue using AMRe87845 at 1:1000.



Immunofluorescence analysis of 3T3 cells labelling LAMP1 with AMRe87845.

注意事项 (Note)

For research use only .