

产品名称: Monoclonal Mouse Anti-Human IgG(Fc specific)
产品货号: AMS0001



产品概述

产品名称 (Production Name)	Monoclonal Mouse Anti-Human IgG(Fc specific)
克隆 (Clonality)	Monoclonal
宿主 (Host)	Mouse
反应种属 (Reactivity)	Human
应用 (Application)	ELISA, WB
形态 (Form)	Liquid
浓度 (Concentration)	5mg/mL
免疫原 (Immunogen)	IgG
特异性 (Specificity)	IgG(Fc specific)
缓冲液 (Buffer)	0.01M PBS、PH:7.2、N preservative
偶联 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

ELISA: 1:20000-1:100000, WB: 1:20000-50000. This monoclonal antibody has 2 clones with good specificity, HG1-3/HG1-11, which only react with human IgG1, IgG2, IgG3, and IgG4, and do not bind to human IgM/IgA/IgE. The heavy chain of HG1-3 antibody is IgM class; the heavy chain of HG1-11 antibody is IgG1 subclass; the light chain is Kappa. Among them, HG1-3 is suitable for coating microplates or nitrocellulose membranes to capture human IgG. The HG1-11 clone is suitable for HRP labeling, etc.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.