

产品名称: Monoclonal Mouse Anti-Human IgM(μ -chain specific)
产品货号: AMS0003



产品概述

产品名称 (Production Name)	Monoclonal Mouse Anti-Human IgM(μ -chain specific)
克隆 (Clonality)	Monoclonal
宿主 (Host)	Mouse
反应种属 (Reactivity)	Human
应用 (Application)	ELISA, WB
形态 (Form)	Liquid
浓度 (Concentration)	5mg/mL
免疫原 (Immunogen)	IgM
特异性 (Specificity)	IgM(μ -chain specific)
缓冲液 (Buffer)	0.01M PBS、PH:7.2、N preservative
偶联 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

ELISA: 1:20000-1:100000, WB: 1:20000-50000. This monoclonal antibody has 2 clones with good specificity, HM-C10/HM-C6, which only react with human IgM and do not bind to human IgG/IgA/IgE. The heavy chains of the antibodies are all IgG1, and the light chains are all Kappa. Among them, HM-C10 is suitable for coating or fixing on agarose microspheres to capture or obtain purified human IgM, and can also be used in related research experiments of sandwich detection of human IgM with HM-C6 labeling.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.