

**产品名称:** BIOTIN-Monoclonal Mouse Anti-Human (Fc specific)  
**产品货号:** AMS0047



## 产品概述

|                        |                                                  |
|------------------------|--------------------------------------------------|
| 产品名称 (Production Name) | BIOTIN-Monoclonal Mouse Anti-Human (Fc specific) |
| 克隆 (Clonality)         | Monoclonal                                       |
| 宿主 (Host)              | Mouse                                            |
| 反应种属 (Reactivity)      | Human                                            |
| 应用 (Application)       | Flow-Cyt/IF/IP,IHC,ELISA,WB                      |
| 形态 (Form)              | Liquid                                           |
| 浓度 (Concentration)     | 1mg/mL                                           |
| 免疫原 (Immunogen)        |                                                  |
| 特异性 (Specificity)      | (Fc specific)                                    |
| 缓冲液 (Buffer)           | 50% Glycerol、0.01M PBS、0.5%BSA、N preservative    |
| 偶联 (Conjugation)       | BIOTIN                                           |
| 修饰 (Modification)      | Unmodification                                   |
| 纯化 (Purification)      | Affinity purification                            |
| 纯度 (Purity)            | Calculated as antibodies, no less than 90%.      |
| 储存方法 (Storage)         | Please store at -20°C away from light.           |

## 使用建议

Flow-Cyt/IF/IP:1:200-2000, IHC: 1:500-2500, ELISA/WB: 1: 2000-20000. This data is a recommended range, please refer to the actual value according to the specific situation.

## 注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.