

产品名称: Monoclonal Mouse Anti-Chicken IgM, μ Chain Specific
产品货号: AMS0503

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	Monoclonal Mouse Anti-Chicken IgM, μ Chain Specific
克隆 (Clonality)	Monoclonal IgG1,Kappa
宿主 (Host)	Mouse
反应种属 (Reactivity)	Chicken
应用 (Application)	ELISA,IF,FC,WB,IHC
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
特异性 (Specificity)	IgM, μ Chain Specific
缓冲液 (Buffer)	0.01M PBS,pH7.2,0.02% Sodium Azide
偶联 (Conjugation)	Unconjugated
最大激发光 (Ex-max)	none
最大发射光 (Em-max)	none
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.

使用建议

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.