

产品名称: Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG)
产品货号: APS0015



产品概述

产品名称 (Production Name)	Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG)
克隆 (Clonality)	Polyclonal
宿主 (Host)	Goat
反应种属 (Reactivity)	Rabbit
应用 (Application)	ELISA, ICC/IF, Flow Cyt, WB
形态 (Form)	Liquid
浓度 (Concentration)	5mg/mL
免疫原 (Immunogen)	IgG
特异性 (Specificity)	IgG (H+L)(Non cross with Human IgG)
缓冲液 (Buffer)	0.01M PBS、PH:7.2、N preservative
偶联 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

ELISA: 1:20000-1:100000; ICC/IF: 1:200-1:2000; Flow Cyt: 1:1000-4000; WB: 1:20000-50000; This data is a recommended range, please refer to the actual value according to the specific situation. This product can also be well applied to conventional in vitro immunological experiments such as colloidal gold.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.