

**产品名称:** HRP-Polyclonal Goat Anti-Mouse IgM( $\mu$ -chain specific)

**产品货号:** APS0020



## 产品概述

|                        |                                                            |
|------------------------|------------------------------------------------------------|
| 产品名称 (Production Name) | HRP-Polyclonal Goat Anti-Mouse IgM( $\mu$ -chain specific) |
| 克隆 (Clonality)         | Polyclonal                                                 |
| 宿主 (Host)              | Goat                                                       |
| 反应种属 (Reactivity)      | Mouse                                                      |
| 应用 (Application)       | ELISA, WB, IHC                                             |
| 形态 (Form)              | Liquid                                                     |
| 浓度 (Concentration)     | 1mg/mL                                                     |
| 免疫原 (Immunogen)        | IgM                                                        |
| 特异性 (Specificity)      | IgM( $\mu$ -chain specific)                                |
| 缓冲液 (Buffer)           | 50% Glycerol、0.01M PBS、Proclin-300                         |
| 偶联 (Conjugation)       | HRP                                                        |
| 修饰 (Modification)      | Unmodification                                             |
| 纯化 (Purification)      | Affinity purification                                      |
| 纯度 (Purity)            | Calculated as antibodies, no less than 90%.                |
| 储存方法 (Storage)         | Please store at -20°C away from light.                     |

## 使用建议

ELISA: 1:2500-5000, WB: 1:2000-4000, This data is a recommended range, please refer to the actual value according to the specific situation. The recommended diluent is 0.01M PBS (PH: 7.0) containing 0.2% BSA and 0.05% Tween-20. Other suitable diluents can also be used. In principle, the prepared working solution should be used up once. If it needs to be stored for a long time and used multiple times after dilution, please contact us to purchase stabilizers or ready-to-use diluents.

## 注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.