

产品名称: HRP-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)
产品货号: APS0029



产品概述

产品名称 (Production Name)	HRP-Polyclonal Goat Anti-Mouse IgG (H+L)(Non cross with Rabbit IgG)
克隆 (Clonality)	Polyclonal
宿主 (Host)	Goat
反应种属 (Reactivity)	Mouse
应用 (Application)	ELISA, WB, IHC
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
免疫原 (Immunogen)	IgG
特异性 (Specificity)	IgG (H+L)(Non cross with Rabbit IgG)
缓冲液 (Buffer)	50% Glycerol、0.01M PBS、Proclin-300
偶联 (Conjugation)	HRP
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

ELISA: 1:4000-20000, WB: 1:4000-12000, IHC: 1:500-2500. This data is a recommended range, please refer to the actual value according to the specific situation. The recommended diluent is 0.01M PBS (PH: 7.0) containing 0.2% BSA and 0.05% Tween-20. Other suitable diluents can also be used. In principle, the prepared working solution should be used up once. If it needs to be stored for a long time and used multiple times after dilution, please contact us to purchase stabilizers or ready-to-use diluents.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.