

产品名称: FITC-Polyclonal Rabbit Anti-Human IgM(μ -chain specific)
产品货号: APS0038



产品概述

产品名称 (Production Name)	FITC-Polyclonal Rabbit Anti-Human IgM(μ -chain specific)
克隆 (Clonality)	Polyclonal
宿主 (Host)	Rabbit
反应种属 (Reactivity)	Human
应用 (Application)	Flow-Cyt/IF,IHC
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
免疫原 (Immunogen)	IgM
特异性 (Specificity)	IgM(μ -chain specific)
缓冲液 (Buffer)	50% Glycerol、0.01M PBS、0.5%BSA、N preservative
偶联 (Conjugation)	FITC
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

Flow-Cyt/IF:1:150-1500, IHC: 1:500-2500, This data is a recommended range, please refer to the actual value according to the specific situation. This product also can be used in Elisa experiment.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.