

产品名称: BIOTIN-Polyclonal Goat Anti-Rabbit IgG (H+L)
(Non cross with Mouse IgG)
产品货号: APS0059



产品概述

产品名称 (Production Name)	BIOTIN-Polyclonal Goat Anti-Rabbit IgG (H+L)(Non cross with Mouse IgG)
克隆 (Clonality)	Polyclonal
宿主 (Host)	Goat
反应种属 (Reactivity)	Rabbit
应用 (Application)	Flow-Cyt/IF/IP,IHC,ELISA,WB
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
免疫原 (Immunogen)	IgG
特异性 (Specificity)	IgG (H+L)(Non cross with Mouse IgG)
缓冲液 (Buffer)	50% Glycerol、0.01M PBS、0.5%GSA、N preservative
偶联 (Conjugation)	BIOTIN
修饰 (Modification)	Unmodification
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.

使用建议

Flow-Cyt/IF/IP:1:200-2000, IHC: 1:500-2500, ELISA/WB: 1: 2000-20000. This data is a recommended range, please refer to the actual value according to the specific situation.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.