

产品名称: R-PE-conjugated Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG)

产品货号: APS0155

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	R-PE-conjugated Polyclonal Goat Anti-Rabbit IgG(H+L)(Non cross with Mouse IgG)
克隆 (Clonality)	Polyclonal Whole IgG
宿主 (Host)	Goat
反应种属 (Reactivity)	Rabbit
应用 (Application)	IF,FC
形态 (Form)	Liquid
浓度 (Concentration)	0.5mg/mL
特异性 (Specificity)	IgG(H+L)(Non cross with Mouse IgG)
缓冲液 (Buffer)	0.01M PBS,0.5%BSA,0.1%PC300,pH7.4
偶联 (Conjugation)	R-PE
最大激发光 (Ex-max)	496nm,565nm
最大发射光 (Em-max)	575nm
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Store at 2-8°C and do not freeze. Prepare working dilution on day of use.Stable for more than six months after shipment.

使用建议

1:100-1:1000.This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.