

产品名称: Polyclonal Goat Anti-Mouse IgM, μ Chain Specific
产品货号: APS0186

Note: Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	Polyclonal Goat Anti-Mouse IgM, μ Chain Specific
克隆 (Clonality)	Polyclonal Whole IgG
宿主 (Host)	Goat
反应种属 (Reactivity)	Mouse
应用 (Application)	ELISA, IF, FC, WB, IHC
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
特异性 (Specificity)	IgM, μ Chain Specific
缓冲液 (Buffer)	0.01M PBS, pH 7.2, 0.02% Sodium Azide
偶联 (Conjugation)	Unconjugated
最大激发光 (Ex-max)	none
最大发射光 (Em-max)	none
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Store at -20°C. Avoid repeated freezing and thawing. Stable for more than one year after shipment.

使用建议

Used in WB, ELISA, IF, IP, IHC, FCM, etc. Before labeling, perform dialysis or ultrafiltration to replace the solution.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.