

产品名称: Fluor750-conjugated Polyclonal Donkey Anti-Goat IgG(H+L) preadsorbed
产品货号: APS0302

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	Fluor750-conjugated Polyclonal Donkey Anti-Goat IgG(H+L) preadsorbed
克隆 (Clonality)	Polyclonal Whole IgG
宿主 (Host)	Donkey
反应种属 (Reactivity)	Goat
应用 (Application)	IF,FC
形态 (Form)	Liquid
浓度 (Concentration)	0.5mg/mL
特异性 (Specificity)	IgG(H+L) Minimal Cross Reactivity:Chicken, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, Rat Serum Proteins
缓冲液 (Buffer)	50% Glycerol,0.01M PBS,0.5%BSA,0.05%PC300,pH7.4
偶联 (Conjugation)	Fluor750
最大激发光 (Ex-max)	747nm
最大发射光 (Em-max)	770nm
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Please store at -20°C away from light.Stable for one year after shipment.

使用建议

1:100-1:500.This data is a recommended range, please refer to the actual value according to the specific situation.
The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.