

产品名称: Polyclonal Donkey Anti-Mouse IgG(H+L) preadsorbed
产品货号: APS0344

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	Polyclonal Donkey Anti-Mouse IgG(H+L) preadsorbed
克隆 (Clonality)	Polyclonal Whole IgG
宿主 (Host)	Donkey
反应种属 (Reactivity)	Mouse
应用 (Application)	ELISA,IF,FC,WB,IHC
形态 (Form)	Liquid
浓度 (Concentration)	1mg/mL
特异性 (Specificity)	IgG(H+L) Minimal Cross Reactivity:Bovine, Chicken, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Rabbit, Sheep Serum Proteins
缓冲液 (Buffer)	0.01M PBS,pH7.2,0.02% Sodium Azide
偶联 (Conjugation)	Unconjugated
最大激发光 (Ex-max)	none
最大发射光 (Em-max)	none
纯化 (Purification)	Affinity purification
纯度 (Purity)	Calculated as antibodies, no less than 90%.
储存方法 (Storage)	Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.

使用建议

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.