

**产品名称:** Polyclonal Donkey Anti-Chicken IgY(IgG)(H+L) preadsorbed  
**产品货号:** APS0380

Note:Centrifuge before opening to ensure complete recovery of vial contents

## 产品概述

|                        |                                                                                                                                         |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| 产品名称 (Production Name) | Polyclonal Donkey Anti-Chicken IgY(IgG)(H+L) preadsorbed                                                                                |
| 克隆 (Clonality)         | Polyclonal Whole IgG                                                                                                                    |
| 宿主 (Host)              | Donkey                                                                                                                                  |
| 反应种属 (Reactivity)      | Chicken                                                                                                                                 |
| 应用 (Application)       | ELISA,IF,FC,WB,IHC                                                                                                                      |
| 形态 (Form)              | Liquid                                                                                                                                  |
| 浓度 (Concentration)     | 1mg/mL                                                                                                                                  |
| 特异性 (Specificity)      | IgY(IgG)(H+L) Minimal Cross Reactivity:Bovine, Goat, Guinea Pig, Syrian Hamster, Horse, Human, Mouse, Rabbit, Rat, Sheep Serum Proteins |
| 缓冲液 (Buffer)           | 0.01M PBS,pH7.2,0.02% Sodium Azide                                                                                                      |
| 偶联 (Conjugation)       | Unconjugated                                                                                                                            |
| 最大激发光 (Ex-max)         | none                                                                                                                                    |
| 最大发射光 (Em-max)         | none                                                                                                                                    |
| 纯化 (Purification)      | Affinity purification                                                                                                                   |
| 纯度 (Purity)            | Calculated as antibodies, no less than 90%.                                                                                             |
| 储存方法 (Storage)         | Store at -20°C.Avoid repeated freezing and thawing. Stable for more than one year after shipment.                                       |

## 使用建议

Used in WB,ELISA,IF,IP,IHC,FCM,etc.Before labeling, perform dialysis or ultrafiltration to replace the solution.

## 注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.