

产品名称: 碱性磷酸酶标记链霉亲和素(SA-ALP)

产品货号: RAD25002

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	碱性磷酸酶标记链霉亲和素(SA-ALP), Alkaline Phosphatase-conjugated Streptavidin
反应种属 (Reactivity)	Biotin
应用 (Application)	IHC/ICC,ELISA,WB
形态 (Form)	Liquid
浓度 (Concentration)	0.5mg/mL
缓冲液 (Buffer)	50% Glycerol,0.005 M Tris-HCl, 0.125 M NaCl, 0.5%BSA,0.02% N preservative,pH 8.0
偶联 (Conjugation)	Alkaline Phosphatase
最大激发光 (Ex-max)	none
最大发射光 (Em-max)	none
储存方法 (Storage)	Please store at -20°C and avoid repeated freezing and thawing .Stable for one year after shipment.

使用建议

IHC/ICC:1:250-500

ELISA、WB:1:500-5000.This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.