

产品名称: PE-Cy5.5标记链霉亲和素

产品货号: RAD25019

Note:Centrifuge before opening to ensure complete recovery of vial contents

产品概述

产品名称 (Production Name)	PE-Cy5.5标记链霉亲和素, PE-Cy5.5-conjugated Streptavidin
反应种属 (Reactivity)	Biotin
应用 (Application)	FC,Immunoassay
形态 (Form)	Liquid
浓度 (Concentration)	0.5mg/mL
缓冲液 (Buffer)	0.01M PBS,0.5%BSA,0.1%Proclin-300,pH 7.4
偶联 (Conjugation)	PE-Cy5.5
最大激发光 (Ex-max)	496nm, 565nm
最大发射光 (Em-max)	696nm
储存方法 (Storage)	Please store at 2-8°C away from light and absolutely avoid freezing.Stable for more than 6 months after shipment.

使用建议

0.1-5 µg/mL (For most applications, Most commonly used in flow cytometry applications) .This data is a recommended range, please refer to the actual value according to the specific situation. The optimal dilution ratio is influenced by various factors, such as antigen concentration and permeability. The actual dilution ratio to be used must be determined through experiments.

注意事项

Avoid freeze / thaw cycles even if this product is very stable. If using a small amount each time, please appropriately divide this antibody into small packages and use one tube each time. This antibody is only for research use, not for diagnosis.