


产品名称: Recombinant Bovine Fetuin A

产品货号: PCB2506

概述 (Summary)

别称 (Alternative Names)	Alpha-2-HS-Glycoprotein; Alpha-2-Z-Globulin; Ba-Alpha-2-Glycoprotein; Fetuin-A; AHSB; FETUA
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤0.5 EU/mg
蛋白构建 (Construction)	Recombinant Bovine Fetuin A is produced by our Mammalian cell expression system and the target gene encoding 119-1359 is expressed. with Tag free
Accession #	P12763
表达宿主 (Host)	Mammalian cell
种属 (Species)	Bovine
预测分子量 (Predicted MW)	36.4 kDa
制剂 (Form)	Lyophilized From PBS,5% mannitol and 0.01% Tween 80, pH7.4
运输方式 (Shipping)	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	12 months at -20°C to -80°C in lyophilized state.6 months at -20°C to -80°C under sterile conditions after reconstitution.7-10 days at 2°C to 8°C under sterile conditions after reconstitution.Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

活性图 (Bioactivity image)



Measured in a cell proliferation assay using B16-F1 cells. Recombinant Bovine Fetuin A stimulates adhesion of B16-F1 cells. The ED50 for this effect is ≤ 100 $\mu\text{g/mL}$.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.