


产品名称: Recombinant Human R-Spondin1

产品货号: PCH2535


概述 (Summary)

别称 (Alternative Names)	RSPO1; R-spondin1; RP11-566C13.1; CRISTIN3; FLJ40906; RSPO Rspo1; R-spondin; Rspodin; RP23-325M14.2; Roof plate-specific spondin-1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤10 EU/mg
蛋白构建 (Construction)	Recombinant Human R-Spondin1 is produced by our Mammalian cell expression system and the target gene encoding Ser21-Ala263 is expressed. with Tag free
Accession #	Q2MKA7
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	26.7 kDa
制剂 (Form)	Lyophilized From PBS,5% mannitol and 0.01% Tween 80, pH7.4
运输方式 (Shipping)	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	12 months at -20°C to -80°C in lyophilized state.6 months at -20°C to -80°C under sterile conditions after reconstitution.7-10 days at 2°C to 8°C under sterile conditions after reconstitution.Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)



活性图 (Bioactivity image)



Measured by its ability to induce Topflash reporter bioactivity in CHO cells in the presence of 50ng/ml wnt3a.
The expected ED50 is <1.0 µg/mL

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.