


产品名称: Recombinant Human TNF- α

产品货号: PCH2540


概述 (Summary)

别称 (Alternative Names)	umor Necrosis Factor; Cachectin; TNF-Alpha; Tumor Necrosis Factor Ligand Superfamily Member 2; TNF-a; TNF; TNFA; TNFSF2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤ 10 EU/mg
蛋白构建 (Construction)	Recombinant Human TNF- α is produced by our Mammalian cell expression system and the target gene encoding Val77-Leu233 is expressed. with Tag free
Accession #	P01375
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	17.4 kDa
制剂 (Form)	Lyophilized From PBS, 5% mannitol and 0.01% Tween 80, pH 7.4
运输方式 (Shipping)	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	12 months at -20°C to -80°C in lyophilized state.6 months at -20°C to -80°C under sterile conditions after reconstitution.7-10 days at 2°C to 8°C under sterile conditions after reconstitution.Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)



活性图 (Bioactivity image)



Measured in a cytotoxicity assay using L-929 mouse fibroblast cells in the presence of the metabolic inhibitor actinomycin D. The ED₅₀ for this effect is 5-60 pg/mL.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.