


产品名称: Recombinant Human bFGF/FGF-2

产品货号: PEH2501


概述 (Summary)

别称 (Alternative Names)	Fibroblast Growth Factor 2; FGF-2; Basic Fibroblast Growth Factor; bFGF; Heparin-Binding Growth Factor 2; HBGF-2; FGF2; FGFB
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤10 EU/mg
蛋白构建 (Construction)	Recombinant Human bFGF/FGF-2 is produced by our E coli expression system and the target gene encoding Pro 143-Ser 288 is expressed. with Tag free
Accession #	P09038
表达宿主 (Host)	E coli
种属 (Species)	Human
预测分子量 (Predicted MW)	17 kDa
制剂 (Form)	Lyophilized From 10mM PB,250mM NaCl,5% mannitol and 0.01% Tween 80,pH 7.4
运输方式 (Shipping)	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	12 months at -20°C to -80°C in lyophilized state.6 months at -20°C to -80°C under sterile conditions after reconstitution.7-10 days at 2°C to 8°C under sterile conditions after reconstitution.Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)



活性图 (Bioactivity image)



Measured in a cell proliferation assay using NIH3T3 mouse embryonic fibroblast Cells. The ED50 for this effect ≤ 2.0 ng/mL.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.