


产品名称: Recombinant Human Insulin

产品货号: PPH2528

概述 (Summary)

别称 (Alternative Names)	Insulin; INS; IDDM; ILPR; IRDN; MODY10
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤10 EU/mg
蛋白构建 (Construction)	Recombinant Human Insulin is produced by our Yeast expression system and the target gene encoding α chain: G-I-V-E-Q-C-C-T-S-I-C-S-L-Y-Q-L-E-N-Y-C-N & beta chain: F-V-N-Q-H-L-C-G-S-H-L-V-E-A-L-Y-L-V-C-G-E-R-G-F-F-Y-T-P-K-T is expressed. with Tag free
Accession #	P01308
表达宿主 (Host)	Yeast
种属 (Species)	Human
预测分子量 (Predicted MW)	6 kDa
制剂 (Form)	Lyophilized From Acetic acid
运输方式 (Shipping)	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	12 months at -20°C to -80°C in lyophilized state.6 months at -20°C to -80°C under sterile conditions after reconstitution.7-10 days at 2°C to 8°C under sterile conditions after reconstitution.Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)



注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.