

产品名称: GMP Recombinant Human Activin A
产品货号: PCH90001




概述 (Summary)


英文全称	Activin A
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤ 10 EU/mg
蛋白构建 (Construction)	Recombinant Human Activin A is produced by our Mammalian cell expression system and the target gene encoding Gly311-Ser426 is expressed.
Accession #	P08476
蛋白标签 (Tag)	Tag free
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	12.9 kDa
蛋白形态 (Form)	Lyophilized
储存缓冲液 (Buffer)	0.085% TFA, 30% ACN, 5% mannitol, pH 2.5
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 8°C under sterile conditions after reconstitution. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: GMP Recombinant Human Activin A
产品货号: PCH90001



活性验证图 (Bioactivity image)



背景 (Background)

分子别名 (Alternative Names)

Activin beta-A chain;Erythroid differentiation protein ;EDF

背景介绍 (References)

Activin and inhibin are two closely related protein complexes that have almost directly opposite biological effects. Activins, members of the TGF-beta superfamily, are disulfide-linked dimeric proteins originally purified from gonadal fluids as proteins that stimulated pituitary follicle stimulating hormone (FSH) release. Inhibins/activins are involved in regulating a number of diverse functions such as hypothalamic and pituitary hormone secretion, gonadal hormone secretion, germ cell development and maturation, erythroid differentiation, insulin secretion, nerve cell survival, embryonic axial development or bone growth, depending on their subunit composition. Activins are homodimers or heterodimers of the various beta subunit isoforms, while inhibins are heterodimers of a unique alpha subunit and one

产品名称: GMP Recombinant Human Activin A
产品货号: PCH90001



of the various beta subunits.

注意事项 (Note)

For research use only .