

产品名称: GMP Recombinant Human Fibronectin
产品货号: PCH90013




概述 (Summary)


英文全称	Fibronectin
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤ 10 EU/mg
蛋白构建 (Construction)	Recombinant Human Fibronectin is produced by our Mammalian cell expression system and the target gene encoding Pro1270-Ser1546&Ala1721-Thr2016 is expressed.
Accession #	P02751
蛋白标签 (Tag)	Tag free
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	62.6 kDa
蛋白形态 (Form)	Lyophilized
储存缓冲液 (Buffer)	PBS, 5% mannitol and 0.01% Tween 80, pH7.4
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 8°C under sterile conditions after reconstitution. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: GMP Recombinant Human Fibronectin
产品货号: PCH90013



活性验证图 (Bioactivity image)



The ED50 for this effect is 35.0-420 ng/mL

背景 (Background)

分子别名 (Alternative Names)

FN; Cold-insoluble globulin; CIG; FN; Fibronectin 1

背景介绍 (References)

Fibronectin1(FN1) is a secreted protein and contains 12 fibronectin type-I domains, fibronectin type-II domains and 16 fibronectin type-III domains. Recombinant human fibronectin fragment, is a protein of ~63 kDa containing a central cell-binding domain, a high affinity heparin-binding domain II, and CS1 site within the alternatively spliced III CS region of human fibronectin. Cells bind to a VLA-4 ligand, a CS-I site, and a VLA-5 ligand, a cell attachment domain, and virus vectors binds to a heparin binding domain II, which co-locates the cell and the virus vector on NovoNectin. This process enhances the density of both cells and vectors, and facilitates the gene transduction in the result.

产品名称: GMP Recombinant Human Fibronectin
产品货号: PCH90013



注意事项 (Note)

For research use only .