

产品名称: GMP Recombinant Human IL-15
产品货号: PCH90022




概述 (Summary)


英文全称	IL-15
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤10 EU/mg
蛋白构建 (Construction)	Recombinant Human IL-15 is produced by our Mammalian cell expression system and the target gene encoding Asn49-Ser162 is expressed.
Accession #	P40933
蛋白标签 (Tag)	Tag free
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	12.7 kDa
蛋白形态 (Form)	Lyophilized
储存缓冲液 (Buffer)	PBS, 5% mannitol and 0.01% Tween 80, pH7.4
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 8°C under sterile conditions after reconstitution. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: GMP Recombinant Human IL-15
产品货号: PCH90022



活性验证图 (Bioactivity image)



The ED50 for this effect is ≤ 0.5 ng/ml

背景 (Background)

分子别名 (Alternative Names)

Interleukin-15; IL-15; IL15

背景介绍 (References)

Human Interleukin 15 (IL-15) is a cytokine that regulates T cell and natural killer cell activation and proliferation. IL-15 binds to the alpha subunit of the IL15 receptor (IL-15RA) with high affinity. IL-15 also binds to the beta and gamma chains of the IL-2 receptor, but not the alpha subunit of the IL2 receptor. IL-15 is structurally and functionally related to IL-2. Both cytokines share some subunits of receptors, allowing them to compete for and negatively regulate each others activity. The number of CD8⁺ memory T cells is controlled by a balance between IL-15 and IL-2. Despite their many overlapping functional properties, IL-2 and IL-15 are, in fact, quite distinct players in the immune system. IL-15 is constitutively expressed by a wide variety of cell types and tissues, including monocytes, macrophages and DCs.

产品名称: GMP Recombinant Human IL-15
产品货号: PCH90022



Mature Human IL-15 shares 70% amino acid sequence identity with Mouse and Rat IL-15.

注意事项 (Note)

For research use only .