

产品名称: GMP Recombinant Human Wnt-3a(hFc Tag)
产品货号: PCH90063




概述 (Summary)


英文全称	Wnt-3a(hFc Tag)
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	≤10 EU/mg
蛋白构建 (Construction)	Recombinant Human Wnt-3a(hFc Tag) is produced by our Mammalian cell expression system and the target gene encoding Ser19-Lys352 with a C-terminal linker GGGGGS and Human IgG1 (Asp104-Lys330) is expressed.
Accession #	NP_149122
蛋白标签 (Tag)	hFc Tag
表达宿主 (Host)	Mammalian cell
种属 (Species)	Human
预测分子量 (Predicted MW)	58.2 kDa
蛋白形态 (Form)	Lyophilized
储存缓冲液 (Buffer)	PBS, pH7.4
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	36 months at -20°C to -80°C in lyophilized state. 6 months at -20°C to -80°C under sterile conditions after reconstitution. 7-10 days at 2°C to 8°C under sterile conditions after reconstitution. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

产品名称: GMP Recombinant Human Wnt-3a(hFc Tag)
产品货号: PCH90063



活性验证图 (Bioactivity image)



The ED50 for this effect is < 30ng/ml..

背景 (Background)

分子别名 (Alternative Names)

MGC119418; MGC119419; MGC119420; protein Wnt-3a; wingless-type MMTV integration site family, member 3A; Wnt3a; Wnt-3a

背景介绍 (References)

Wnt surrogate is a Wnt pathway agonists and related compositions, which may be used in any of a variety of therapeutic methods for the treatment of diseases.

注意事项 (Note)

For research use only .