


Product Name: Recombinant Human AKR1C2
Catalog #: PEH0039

概述 (Summary)

英文全称	AKR1C2/Aldo-keto reductase family 1 member C2
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Aldo-Keto Reductase Family 1 Member C2 is produced by our E.coli expression system and the target gene encoding Met1-Tyr323 is expressed.
Accession #	P52895
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	36.74 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human AKR1C2
Catalog #: PEH0039



分子别名 (Alternative Names)

Aldo-Keto Reductase Family 1 Member C2; 3-Alpha-HSD3; Chlordecone Reductase Homolog HAKRD; Dihydrodiol Dehydrogenase 2; DD-2; DD2; Dihydrodiol Dehydrogenase/Bile Acid-Binding Protein; DD/BABP; Trans-1; 2-Dihydrobenzene-1; 2-Diol Dehydrogenase; Type III 3-Alpha-Hydroxysteroid Dehydrogenase; AKR1C2; DDH2

背景介绍 (References)

Aldo-Keto Reductase Family 1 Member C2 (AKR1C2) plays a role in concert with the 5- α /5- β -Steroid Reductases to convert Steroid hormones into the 3- α /5- α and 3- α /5- β -Tetrahydrosteroids. AKR1C2 catalyzes the inactivation of the most potent androgen 5- α -Dihydrotestosterone (5- α -DHT) to 5- α -Androstane-3- α , 17- β -diol (3- α -diol).

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.