


Product Name: Recombinant Human GALM (C-6His)
Catalog #: PEH0040



概述 (Summary)

英文全称	Aldose 1-epimerase/GALM
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Galactose Mutarotase is produced by our E.coli expression system and the target gene encoding Ala2-Ala342 is expressed with a 6His tag at the C-terminus.
Accession #	Q96C23
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	38.8 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names)	Aldose 1-Epimerase; Galactose Mutarotase; GALM
--------------------------	--

Product Name: Recombinant Human GALM (C-6His)
Catalog #: PEH0040



背景介绍 (References)

Galactose Mutarotase (GALM) is a cytoplasmic enzyme that belongs to the Aldose Epimerase family. GALM is a Mutarotase that converts α -Aldose to the β -Anomer. GALM is active on D-Glucose, L-Arabinose, D-Xylose, D-Galactose, Maltose and Lactose. GALM may be required for normal Galactose metabolism by maintaining the equilibrium of α - and β - anomers of Galactose.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.