


Product Name: Recombinant Human ASB13 (N-6His)
Catalog #: PEH0067



概述 (Summary)

英文全称	Ankyrin repeat and SOCS box protein 13/ASB13
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ankyrin Repeat And SOCS Box Protein 13 is produced by our E.coli expression system and the target gene encoding Met1-Asn278 is expressed with a 6His tag at the N-terminus.
Accession #	Q8W XK3
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	31.4 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, 4M Urea, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Ankyrin repeat and SOCS box protein 13; ASB13;

Product Name: Recombinant Human ASB13 (N-6His)
Catalog #: PEH0067



背景介绍 (References)

Ankyrin repeat and SOCS box protein 13(ASB13) is a member of the ankyrin repeat and SOCS box-containing (ASB) family. ASB13 contain six ankyrin repeats sequence and a SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation. ASB13 may be a substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.