


概述 (Summary)

英文全称	Annexin A8/ANXA8
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Annexin A8 is produced by our E.coli expression system and the target gene encoding Met1-Pro327 is expressed.
Accession #	P13928
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	36.9 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Lyophilized protein should be stored at ≤ -20°C, stable for one year after receipt. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at ≤ -20°C for 3 months.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human ANXA8
Catalog #: PEH0076



背景 (Background)

分子别名 (Alternative Names)

Annexin A8; Annexin VIII; Annexin-8; Vascular Anticoagulant-Beta; VAC-Beta; ANXA8; ANX8

背景介绍 (References)

Annexin A8 (ANXA8) belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Annexin A8 contains 4 annexin repeats separated by linking sequences of variable lengths. A pair of annexin repeats may form one binding site for calcium and phospholipid. ANXA8 is preferentially expressed in acute promyelocytic leukemia (APL) cells, which suggests its participation in hematopoietic cell differentiation. In addition, ANXA8 is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.