


Product Name: Recombinant Human ATG3
Catalog #: PEH0110

概述 (Summary)

英文全称	ATG3/Apg3L/Apg3p
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Ubiquitin-like-conjugating Enzyme ATG3/ATG3 is produced by our E.coli expression system and the target gene encoding Met1-Met314 is expressed.
Accession #	Q9NT62
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	36 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 150mM NaCl, 10% Glycerol, pH 8.0.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human ATG3
Catalog #: PEH0110

分子别名 (Alternative Names)

Ubiquitin-like-conjugating enzyme ATG3; Autophagy-related protein 3; APG3-like; hApg3; Protein PC3-96; ATG3; APG3; APG3L

背景介绍 (References)

Ubiquitin-like-conjugating enzyme ATG3 (ATG3), also known as Apg3L and Apg3p, functions as a regulatory component of autophagosome biogenesis necessary for autophagy. ATG3 exhibits 98% aa sequence identity with both its mouse and rat orthologs. It is widely expressed and has highly levels in heart, skeletal muscle, kidney, liver and placenta. As an E2-like enzyme, involves in autophagy and mitochondrial homeostasis. ATG3 catalyzes the conjugation of ATG8-like proteins to PE which is essential for autophagy. ATG3 also can catalyze the conjugation of ATG12 to itself which plays a role in mitochondrial homeostasis but not in autophagy.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.