


产品名称: Recombinant Human CASP10 (C-6His)
产品货号: PEH0236



概述 (Summary)

英文全称	Caspase-10/CASP10
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Caspase-10 is produced by our E.coli expression system and the target gene encoding Val220-Ile480 is expressed with a 6His tag at the C-terminus.
Accession #	Q92851-4
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	30.1 KDa
制剂 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 8% Sucrose, 1mM DTT, 0.05% Tween80, pH8.5.
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

产品名称: Recombinant Human CASP10 (C-6His)
产品货号: PEH0236



分子别名 (Alternative Names)

Caspase-10; CASP-10; Apoptotic Protease Mch-4; FAS-Associated Death Domain Protein Interleukin-1B-Converting Enzyme 2; FLICE2; ICE-Like Apoptotic Protease 4; CASP10; MCH4

背景介绍 (References)

Caspase-10 (CASP10) is a 521 amino acid protein member of the Cysteine-Aspartic Acid Protease (Caspase) family. CASP10 contains two DED (Death Effector) domains and is detectable in most tissues. CASP10 cleavage by Granzyme B and autocatalytic activity generate the two active subunits: Caspase-10 subunit p23/17, Caspase-10 subunit p12. Caspases are a family of cytosolic aspartate-specific cysteine proteases involved in the execution-phase of cell apoptosis, the initiation and execution. Human caspases can be subdivided into three functional groups: cytokine activation (caspase-1, -4, -5, and -13), apoptosis initiation (caspase-2, -8, -9, -and -10), and apoptosis execution (caspase-3, -6, and -7). CASP10 cleaves and activates caspases 3 and 7, but itself is processed by caspase 8. Defects in CASP10 are associated with apoptosis defects seen in type II autoimmune lymphoproliferative syndrome.

注意事项 (Note)

For research use only .