


Product Name: Recombinant Human CLIC3 (C-6His)
Catalog #: PEH0417



概述 (Summary)

英文全称	CLIC3/Chloride intracellular channel protein 3
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Chloride Intracellular Channel Protein 3 is produced by our E.coli expression system and the target gene encoding Met1-Arg236 is expressed with a 6His tag at the C-terminus.
Accession #	O95833
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	27.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of 20mM Tris-HCl, 200mM Sucrose, 0.05% Tween 80 (w/v), pH 8.5.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

Product Name: Recombinant Human CLIC3 (C-6His)
Catalog #: PEH0417



分子别名 (Alternative Names)

Chloride intracellular channel protein 3; CLIC3

背景介绍 (References)

Chloride intracellular channel protein 3 (CLIC3) is encoded by the CLIC3 gene. CLIC3 is a single-pass membrane protein which belongs to the chloride channel CLIC family. It contains one GST C-terminal domain and one GST N-terminal domain. Chloride intracellular channel protein 3 high expressed in the placental, lung and heart, low expressed in skeletal muscle, kidney and pancreas. Chloride intracellular channel protein 3 can insert into membranes and forms chloride ion channels, may participate in cellular growth control.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.