

Product Name: Recombinant Human DCUN1D1 (N-6His)
Catalog #: PEH0520




概述 (Summary)

英文全称	DCN1-like protein 1/DCUN1D1
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human DCN1-Like Protein 1 is produced by our E.coli expression system and the target gene encoding Met1-Val259 is expressed with a 6His tag at the N-terminus.
Accession #	Q96GG9
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	32.3 KDa
蛋白形态 (Form)	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles. Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100μg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

电泳图 (SDS-PAGE image)

Product Name: Recombinant Human DCUN1D1 (N-6His)
Catalog #: PEH0520



背景 (Background)

分子别名 (Alternative Names)

DCN1-Like Protein 1; DCUN1 Domain-Containing Protein 1; Defective in Cullin Neddylation Protein 1-Like Protein 1; Squamous Cell Carcinoma-Related Oncogene; DCUN1D1; DCUN1L1; RP42; SCCRO

背景介绍 (References)

DCN1-Like Protein 1 is a protein containing 1 DCUN1 domain and 1 UBA-like domain. DCN1-Like Protein 1 contains part of an E3 Ubiquitin Ligase Complex for Neddylation. It is required for Neddylation of Cullin components of E3 Cullin-RING Ubiquitin Ligase complexes. It enhances the rate of Cullins Neddylation. DCN1-Like Protein 1 recruits the NEDD8-charged E2 Enzyme to the Cullin component. DCUN1D1 is involved in the release of inhibitory effects of CAND1 on the Cullin-RING Ligase E3 complex assembly and activity. It acts also as an oncogene facilitating malignant transformation and carcinogenic progression.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.