


Product Name: Recombinant Human dUTPase
Catalog #: PEH0546

概述 (Summary)

英文全称	dUTPase/DUT
纯度 (Purity)	Greater than 95% as determined by reducing SDS-PAGE
内毒素 (Endotoxin level)	<1 EU/μg as determined by LAL test.
蛋白构建 (Construction)	Recombinant Human Deoxyuridine 5-Triphosphate Nucleotidohydrolase, Mitochondrial is produced by our E.coli expression system and the target gene encoding Met1-Asn164 is expressed.
Accession #	P33316-2
蛋白标签 (Tag)	
表达宿主 (Host)	E.coli
种属 (Species)	Human
预测分子量 (Predicted MW)	17.7 KDa
蛋白形态 (Form)	Supplied as a 0.2 μm filtered solution of PBS, pH 7.4.
储存缓冲液 (Buffer)	
运输方式 (Shipping)	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
稳定性&储存 (Stability &Storage)	Store at ≤-70°C, stable for 6 months after receipt. Store at ≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
复溶 (Reconstitution)	

电泳图 (SDS-PAGE image)



背景 (Background)

分子别名 (Alternative Names) Deoxyuridine 5-Triphosphate Nucleotidohydrolase Mitochondrial; dUTPase;

Product Name: Recombinant Human dUTPase
Catalog #: PEH0546



背景介绍 (References)

dUTP Pyrophosphatase; DUT

Deoxyuridine 5-Triphosphate Nucleotidohydrolase Mitochondrial (dUTPase) belongs to the dUTPase family. dUTPase exists as a homotrimer and is involved in nucleotide metabolism. dUTPase produces dUMP, the immediate precursor of thymidine nucleotides and it decreases the intracellular concentration of dUTP so that uracil cannot be incorporated into DNA. The dUTPase increase in PCR product yield, length and fidelity enables further down-stream applications. These effects make dUTPase useful in PCR fidelity and yield-sensitive applications. dUTPase is specific for dUTP and is critical for the fidelity of DNA replication and repair.

注意事项 (Note)

For Research Use Only , Not for Diagnostic Use.